

ONE FAMILY OF
ANALYZERS.

ONE STANDARD
OF TESTING.



One standardized analyzer line. Many hemostasis testing needs.

The **ACL TOP Family** is a complete line of Hemostasis Testing Systems—each designed with the same high standards of operation and functionality to meet the needs of your hospital and lab regardless of its size or scope. Whether running routine or highly specialized assays, ACL TOP provides **Testing Process Automation** in a range of instruments, ideal for low- to high-volume labs through **fully integrated Laboratory Automation Systems**.

And, because ACL TOP Systems are optimized to run with HemosIL reagents, they are part of a **complete disease state management solution**.

All ACL TOP Systems provide:

- **Same results**—standardization among labs
- **Same reagents and consumables**—optimized administration
- **Same features and usability**—simplified training
- **Same powerful and intuitive software**—easy to use

For low- to medium-volume labs

ACL TOP 300 CTS*

Newest and smallest ACL TOP analyzer for routine or specialty assays in low- to medium-volume labs. Perfect companion to ACL TOP 700 CTS or ACL TOP 500 CTS.

Throughput (PT/h)	110
Samples onboard	40
Reagents onboard	26



New

For medium- to high-volume labs

ACL TOP 500 CTS

Highly automated testing processes in routine, medium- to high-volume and specialty labs.

Throughput (PT/h)	240
Samples onboard	80
Reagents onboard	40



For high-volume labs

ACL TOP 700

Run more tests in less time—user-friendly and high throughput for routine analysis in laboratories with the heaviest workloads.

ACL TOP 700 CTS

Added safety in routine, high-volume and specialty labs with Closed-Tube Sampling.

ACL TOP 700 LAS

Maximum flexibility and efficiency; connects to your laboratory automation track.

	ACL TOP 700 CTS	ACL TOP 700
Throughput (PT/h)	270	360
Samples onboard	120	120
Reagents onboard	60	60



*Not available in all countries.

'Testing Process Automation' and superior performance across all ACL TOP Systems

Process more samples in less time with minimal operator intervention.

Analyzer Automation

- Barcoded reagents
- Continuous onboard reagent stability monitoring
- Automatic QC frequency execution
- Rerun and reflex testing capabilities
- Fully automated reporting of factor assays with Parallelism
- Auto-verification and uploading of results
- Closed-Tube Sampling via cap-piercing (CTS models only)

Laboratory Automation

- Follows CLSI Guidelines (AUTO 1–5) for true 'Point-of-Reference' sampling
- Open system compatible with most Laboratory Automation Track Systems
- Eliminates need for costly and slow robotic interface

Continuous Operation

- Continual loading and unloading of samples and reagents via racks with no system interruption
- Uninterrupted cuvette loading
- Uninterrupted waste disposal

Fast Turnaround

- Up to 360 PT/hour
- PT from standby in < 3 minutes
- Samples loaded on any rack, in any position, at any time, including STAT

Simple Maintenance

- Daily maintenance, ordered by user and performed by system, in < 5 minutes
- 'Maintenance overdue' notifications alert user
- Remote instrument diagnostics and troubleshooting via web interface in real time* (optional)

*Not available in all countries.

ACL TOP Systems provide a variety of control and visibility features for enhanced performance and quality reporting—ensuring ultimate testing confidence.

Clot Signature Curves

- Available for patient samples, QC and calibrations
- Illustrate reaction readings within a cuvette, for analysis of abnormal results
- Unique to optical-based detection systems

State-of-the-Art Factor Assay Testing

- Multi-point calibration with direct dilution
- One calibration for entire linear range
- Dual calibration algorithms, including polynomial math models
- Zero calibration point for highest accuracy at very low factor activity
- Automatic multiple dilutions for sample Factor Parallelism
- Graphical display and data checks enhance identification of factor inhibitors



User-friendly interface with large touchscreen is simple to operate.



Closed-Tube Sampling ensures safety.




Onboard barcode reader allows real-time reagent management.



ACL TOP 700 LAS easily integrates into most laboratory automation track systems.

Advanced optical system offers superior accuracy and efficiency

- Wavelengths optimized for minimal interferences from known contaminants
- Multiple- and independent-reading channels
- Simultaneous reading of multiple wavelengths
- Clot curve outputs for independent verification and validation of results



ACL TOP Systems read at 671 nm, outside the range of optical interferences for several key and routine assays.

HemosIL® Reagents

Comprehensive automated test panel for complete disease state management

ACL TOP Family Systems offer complete disease-state management through a comprehensive panel of HemosIL assays. The broad test menu provides flexible solutions for both routine and specialty testing. The same reagents on all ACL TOP Systems mean true standardization throughout the lab.

General Screening and Anticoagulant Monitoring

PT
 APTT
 Fibrinogen Clauss
 Thrombin Time
 Hepatocomplex*
 Pro-IL-Complex*
 Liquid Heparin (one-stage)
 Heparin (two-stage)

D-Dimer

D-Dimer HS 500
 D-Dimer HS
 D-Dimer 500*
 D-Dimer

Heparin-Induced Thrombocytopenia

HIT-Ab_(PF4-H)*

Antiphospholipid Syndrome

Silica Clotting Time
 LAC Screen/Confirm
 dRVVT Screen/Confirm

Coagulation Factors

Intrinsic Factors
 Extrinsic Factors
 FXIII Antigen

von Willebrand Disease

VWF Antigen
 VWF Activity
 VWF:RCo*

Thrombophilia

Antithrombin
 Protein C (Chromogenic)
 Protein C (Clotting)
 Free Protein S (Antigenic)
 Free Protein S (Clotting)
 FV Leiden (APC-R V)
 Homocysteine
 ThromboPath*

Fibrinolysis

Plasminogen
 Plasmin Inhibitor

*Not currently 510(k) cleared.

HemosIL®





Fully customizable for labs of any size

- User-configurable settings
- Cap-piercing option
- Laboratory Automation System option

At IL, we are passionate about bringing the most innovative solutions to hemostasis testing. Through our ACL products, we offer a variety of the highest-quality instruments, all complemented by a full panel of HemosIL assays and supported by IL's worldwide technical team. But more than just a broad hemostasis testing offering, we are committed to continuous enhancements through constant innovation.

▶ HEMOSTASIS INNOVATION IS HERE.

ACL TOP Family Technical Specifications

Test menu		Clotting, chromogenic and immunological assays
System configurations	ACL TOP 700 ACL TOP 700 CTS ACL TOP 700 LAS ACL TOP 500 CTS ACL TOP 300 CTS	Open-Tube Sampling Closed-Tube Sampling (cap-piercing) Laboratory Automation System compatibility (off-board sampling) Compact size with Closed-Tube Sampling Compact size with Closed-Tube Sampling
Samples onboard	ACL TOP 700 ACL TOP 700 CTS ACL TOP 700 LAS ACL TOP 500 CTS ACL TOP 300 CTS	120 120 Continuous from LAS track or 90 front-loaded 80 40
Reagents onboard	ACL TOP 700 ACL TOP 700 CTS ACL TOP 700 LAS ACL TOP 500 CTS ACL TOP 300 CTS	60 60 60 40 26
Patient samples results database		20,000 samples
Throughput	ACL TOP 700 ACL TOP 700 CTS ACL TOP 700 LAS ACL TOP 500 CTS ACL TOP 300 CTS	PT..... 360 tests/hr (360 samples/hr) APTT..... 320 tests/hr (320 samples/hr) PT and APTT 330 tests/hr (165 samples/hr) PT..... 270 tests/hr (270 samples/hr) APTT..... 270 tests/hr (270 samples/hr) PT and APTT 260 tests/hr (130 samples/hr) Throughput from the LAS Track is dependent upon factors such as LAS Track configuration and settings PT..... 240 tests/hr (240 samples/hr) APTT..... 180 tests/hr (180 samples/hr) PT and APTT 180 tests/hr (90 samples/hr) PT..... 110 tests/hr (110 samples/hr) APTT..... 110 tests/hr (110 samples/hr) PT and APTT 110 tests/hr (55 samples/hr)
Dimensions (w x d x h)	ACL TOP 700 ACL TOP 700 CTS ACL TOP 700 LAS ACL TOP 500 CTS ACL TOP 300 CTS	151 x 76 x 73 cm 151 x 76 x 73 cm 188 x 87 x 162 cm (includes integrated table) 110 x 82 x 73 cm 81 x 84 x 73 cm
Weight	ACL TOP 700 ACL TOP 700 CTS ACL TOP 700 LAS ACL TOP 500 CTS ACL TOP 300 CTS	150 kg 160 kg 184 kg (includes integrated table) 141 kg 91 kg

Werfen Group

Worldwide Locations

Corporate Headquarters

Barcelona, Spain
Tel. +34-93-4010101
www.werfengroup.com

Instrumentation Laboratory Headquarters

Bedford, MA
Tel. +1-781-861-0710
www.ilwww.com

US, Canada, Latin America and South America

IL USA
Bedford, MA
Tel. +1-781-861-0710
www.ilus.com

Werfen Medical Ltd. Brazil
São Paulo
Tel. +55-11-4154-3337

IL Canada
Richmond Hill, ON
Tel. +1-800-552-2025
x6115

IL Mexico
Col. Granada
Tel. +52-55-5262-1760
www.il-mexico.com.mx

Izasa Uruguay
Montevideo
Tel. +59-82-4818133

Pacific

Werfen Medical China
Shanghai
Tel. +86-21-32100745

Beijing
Tel. +86-10-359006230

Werfen Hong Kong
Hong Kong
Tel. +852-27927773

IL India
New Delhi
Tel. +91-11-25191251

IL Japan
Minato-ku, Tokyo
Tel. +81-3-3437-6350

Werfen Medical Korea
Seoul
Tel. +82-2-571-9246
www.werfenmedical.com

Europe, Middle East, Africa

IL Austria
Vienna
Tel. +43-1-2565800-0

IL Belgium
Zaventem
Tel. +32-2-7252052
www.il-be.com

Comesa Czech
Prague
Tel. +420-2-7816047

IL France
Le Pré Saint Gervais
Tel. +33(0)1-82-308600
www.il-france.fr

IL Germany
Munich
Tel. +49-89-909070
www.il-ger.de

IL The Netherlands
Breda
Tel. +31(0)-76-5480100
www.il-nl.com

Comesa Hungary
Budapest
Tel. +36-1-4392910
or 11

IL Italy
Milan
Tel. +39-02-25221
www.il-italia.it

IL Lithuania
Kaunas
Tel. +370-37-313157

Comesa Poland
Warsaw
Tel. +48-22-3361800

Izasa Portugal
Carnaxide
Tel. +351-21-4247300
www.izasa.com

IL Russia
Moscow
Tel. +7-495-1244559

Izasa Spain
Barcelona
Tel. +34-93-4010101
www.izasa.com

IL UK
Warrington, Cheshire
Tel. +44-1925-81-0141
www.il-uk.com

Instrumentation Laboratory, Comesa and Izasa are companies of Werfen Group.

For more information, contact your local IL sales representative or distributor.

The Instrumentation Laboratory logo, ACL, ACL TOP and HemosIL are trademarks of Instrumentation Laboratory Company and/or one of its subsidiaries or parent companies, and may be registered in the United States Patent and Trademark Office and in other jurisdictions. All other product names, company names, marks, logos and symbols are trademarks of their respective owners.

©2012 Instrumentation Laboratory. All rights reserved.

