


ELECTROPHORESIS AUTOMATED SYSTEM

ADALYA 40




 **SELEO**
Engineering

Seleo Engineering s.r.l.
1ª traversa bugnano s.c.
81030 Orta di Atella (CE) Italy
TEL +39 (0)81 5022877 FAX +39(0)81 502287702
Homepage www.seleo.it
e-mail info@seleo.it



ELECTROPHORESIS AUTOMATED SYSTEM



FAST/VERSATILE

4/40 samples, up to 5
electrophoresis membranes
8 samples each.

Allows to perform the most common
techniques:

- Serum Proteins
- Haemoglobins
- Lipoproteins
- Urinary Proteins

COMPLETE AUTOMATION

Adalya automatically arranges
the loading of all the reagents
to standardize the delicate
electrophoretic process.

During test all the liquid are
constantly measured and
adjusted, minimising the
variables introduced by reagents
level variations and
to ensure the test preparation in true
complete automation.



Serum Proteins



Haemoglobins



Lipoproteins

PRODUCTIVITY

Serum-Proteins : 24 samples/ 1h:15'

Lipoproteins: 16 samples/ 1h:15'

Haemoglobins: 16 samples/ 1h:15'


SOFTWARE FEATURES

User interface: monitor, keyboard and mouse

Minima recognition: automatic

Managing techniques: adjustable techniques in times, migration voltage/current

Operating system:



ADALYA 40



SELEO
Engineering

Your partner in Protein Electrophoresis