

ELECTROPHORESIS READER

 **MICROPHOR08**



Pocket Electrophoresis Densitometer

- Ready for transparent/opaque scanning
- Fast and reliable
- Suitable for agarose gel and cellulose acetate

 **SELEO**
Engineering

Seleo Engineering s.r.l.
1ª traversa bugnano s.c.
81030 Orta di Atella (CE) Italy
TEL +39 (0)81 5022877 FAX +39(0)81 502287702
Homepage www.seleo.it
e-mail info@seleo.it

**MADE IN
ITALY**

ELECTROPHORESIS READER

Smart and reliable



Technical features

- LED light source and high gain phototransistors
- Reading system in absorbance with self-regulation of the source bright
- Reading media in cellulose acetate (transparent and opaque)
- Reading media in agarose gel
- 8 Reading samples at a time
- Connection to pc via USB
- External power supply 24v
- Dimension: 11.6x15.6x6.3 cm (LxPxH)

Software features

- Pre-programmed Analysis: Serum proteins, Hemoglobins, Lipoproteins
 - Possibility to define additional exams
- Automatic Processing of electrophoretic pattern
 - Possibility of correction of paths
- Personal management and printing reports
 - Report header customization

Operating system :



MICROPHOR08



SELEO
Engineering



Your partner in Protein Electrophoresis